



Beige wrist restraint h983_14

Merino wool wrist restraint with cotton straps. Strap length: 140 cm.

The **H983_14 wrist restraint** is a non-sterile medical device, designed to **safely immobilize the patient's wrist**, for example to prevent the self-removal of devices such as IV drips or catheters. It is indicated for use in **hospitals, nursing homes, intensive care units and home care**, where it is necessary to temporarily limit the patient's movements, always under medical advice.

The wrist restraint is made of **100% Merino wool**, a soft and natural material, ideal for **direct contact with the skin**, even for long periods, as it ensures **comfort, breathability and absence of irritation**. The **straps (braid tapes)** are made of **100% cotton, 140 cm** long, designed for secure attachment to the bed frame or to two different anchor points.

The design of the device allows for **simple and stable fixation**: the tie passes through a metal ring, is knotted securely and allows for optimal movement control. The wrist restraint can be used alone or in combination with other restraint aids. After use, it can be **machine washed at 30°C** with a wool program, avoiding ironing, tumble drying and dry cleaning.

Buy the **H983_14 wrist restraint** now to ensure safety, comfort and control during hospital or home care. Contact us for healthcare supplies in bulk.

Technical specifications:

- **Usage:** wrist containment in a secure position
- **Padding material:** 100% Merino wool
- **Straps (braid tapes):** 100% cotton
- **Strap length:** 140 cm
- **Color:** beige
- **Device:** non-sterile
- **Indicated:** for direct skin contact
- **Washing:** in lukewarm water max 30°C – wool program

User benefits:

- **Prevents:** unwanted actions such as the removal of IV drips or catheters.
- **Material:** soft and natural, safe in contact with the skin.
- **Easy:** to apply and to fix to the bed frame.
- **Straps:** long for multiple or double-point fixation.
- **Washable:** reusable for multiple use cycles.

FAQ:

1. **What is the H983_14 wrist restraint used for?**
It is used to safely restrain the patient's wrist and prevent involuntary movements.
2. **Is it safe to use on bare skin?**
Yes, it is made of natural Merino wool and can be applied directly to the skin.
3. **Are the straps long enough for double fixation?**
Yes, they measure 140 cm and can be fixed to two points of the frame.

4. **How is the device washed?**

It can be machine washed at 30°C, using a delicate or wool program. Do not iron or tumble dry.

5. **Is it suitable for all patients?**

Yes, unless otherwise contraindicated by medical staff. Always follow medical instructions.

Quick buying guide:

When choosing a wrist restraint for medical use, ensure it is **made with hypoallergenic and breathable materials**, such as **Merino wool**.

The **length of the straps** is essential to guarantee effective fixation even on beds of different sizes.

The H983_14 wrist restraint is a reliable solution for clinical or home settings where **temporary patient containment** is necessary.

It must always be used under **medical supervision** and in **compliance with current health regulations**.

Hospital wrist restraint, patient containment wristband, merino wool wrist holder, wrist containment strap, wrist immobilization aid, cotton wrist restraint, IV wrist holder, medical strap for hospital bed.

**Image for illustrative purposes only.*

INFORMATION

