



Electric elliptical pedal exerciser for rehabilitation h1017_06

Mini electric elliptical for leg rehabilitation with remote control and display. Dimensions: approx. W 50 x D 44 x H 30 cm.

The **h1017_06 electric elliptical pedal exerciser** is a compact and effective device, designed to perform **passive and active leg rehabilitation exercises** while **seated**. Ideal for **people with reduced mobility**, it promotes the recovery of motor functions, improves **circulation**, stimulates balance, and helps **strengthen lower limb muscles** in total safety.

This pedal exerciser is equipped with a **multi-function remote control** (requires 2 AAA batteries, not included), which allows you to select **modes, preset programs, movement direction, speed**, and to turn the device on/off. The integrated **LED display** shows real-time **time, speed, distance, calories burned, and rotations**.

The device is mains-powered and offers a **power of 50W**, with an **adjustable speed range from 1 to 30 levels** to adapt to every need. The **two rotation directions (forward and backward)** allow for a complete workout. The **automatic modes P1, P2, and P3**, each with different durations and sequences of speed and direction, make the training varied and progressive.

Correct use of the device, while **seated on a stable chair or armchair**, allows for real physical and cardiovascular benefits. The exercise **must never be performed standing up**, and it is essential to follow the **safety instructions** to ensure proper and safe use over time.

Choose the **h1017_06 electric elliptical pedal exerciser** and improve your mobility safely, directly from home. Contact us for offers on quantities or professional use.

Technical features:

- **Type:** seated electric elliptical pedal exerciser
- **Intended use:** leg exercise and rehabilitation
- **Instructions for use:** seated only (do not stand on the device)
- **LED Display:** shows speed, time, distance, calories, rotations
- **Remote control:** included (requires 2 AAA batteries, not included)
- **Movement directions:** forward and backward
- **Adjustable speed:** 1–30 levels
- **Automatic programs:** P1 (30 min), P2 (15 min), P3 (30 min)
- **Voltage:** 110–220 V
- **Power:** 50 W
- **Frequency:** 50–60 Hz
- **Net weight:** 6.7 kg
- **Dimensions:** approx. W 50 x D 44 x H 30 cm
- **Mains plug:** rear connection, cable included
- **Transport handle:** integrated

User benefits:

- **Promotes:** motor rehabilitation of the legs in a safe way.
- **Improves:** blood circulation and cardiovascular activity.
- **Offers:** passive exercise while seated, also ideal for the elderly.
- **Display:** intuitive to keep track of workout data.

- **Easy:** to use thanks to the remote control with all functions.

FAQ:

1. **Can it be used standing up?**

No. It is mandatory to use the device while seated. Never stand on the pedal exerciser.

2. **How much does the product weigh?**

The product weighs 6.7 kg.

3. **Is it suitable for people with mobility difficulties?**

Yes, it is designed specifically for users with reduced mobility or in the rehabilitation phase.

4. **How does the remote control work?**

The remote control allows you to start/stop the activity, choose direction, mode, speed, and programs.

5. **What are the functions of the display?**

The display shows: time, distance, speed, calories, number of rotations, mode, and active program.

Mini buying guide:

The h1017_06 pedal exerciser is designed for those who need to perform **light and safe physical activity** while seated.

Before purchasing, evaluate **the intended use** (rehabilitation, motor maintenance, prevention of a sedentary lifestyle) and the available space.

Thanks to the **adjustable speed**, the presence of **automatic programs**, and the **remote control**, it is ideal for users who desire **autonomy and simplicity of use**.

It is essential **not to use it standing up**, to fix it on a **flat surface free of obstacles**, and to always disconnect it after use for **greater safety**.

Rehabilitation elliptical pedal exerciser, mini electric elliptical, seated pedal exerciser for the elderly, compact electric elliptical, leg motor rehabilitation, motorized leg pedal exerciser, passive rehabilitation tool.

**Images are purely indicative.*

INFORMATION



HOLITY.COM



HOLITY.COM



HOLITY.COM