



Pack of 2 pieces bed fixation belt h983_03

Spreader strap for leg containment in a spread position, with anti-friction padding and velcro closure. Total length: 260 cm. Adjustable circumference: from 15 to 26 cm. To be used exclusively by healthcare personnel.

The **bed fixation belt H983_03** is a **spreader strap** designed to hold **the ankles in a spread position**, preventing the muscles from working at full power. This soft restraint system is ideal for use in **hospitals, clinics, and healthcare facilities**, ensuring **safe containment** while respecting the patient's comfort.

Made of **100% cotton**, it is equipped with **anti-friction padding** that reduces the risk of skin irritation, while the **velcro closure** allows for practical and fast adjustment.

The total length is **260 cm**, with **six adjustment eyelets in Nickel/Brass alloy on each side**, ensuring **wide adaptability to the patient's needs**.

To be used exclusively by **specialized healthcare personnel**.

For an **additional level of safety**, it is possible to purchase **the safety buttons and the safety key separately (see optional)**, which prevent unauthorized opening of the device. Removal is possible **only via the dedicated magnetic key**.

Buy the **H983_03 spreader strap** now for safe and adjustable containment in healthcare contexts! Request a personalized quote immediately.

Technical features:

- **Material:** 100% cotton with anti-friction padding
- **Closure:** Adjustable velcro
- **Adjustment:** Six Nickel/Brass alloy eyelets per side
- **Total length:** 260 cm
- **Adjustable circumference:** from 15 to 26 cm
- **Color:** Cream
- **Use:** Limitation of leg movement in a spread position
- **Safety option:** Safety buttons and safety key purchasable separately (see optional)
- **Pack:** of two pieces

Optional:

- **2 safety buttons:** automatic closing buttons to be inserted into the overlapping eyelets of the safety belts
- **1 safety key:** polyamide key with reinforced neodymium magnet 8 mm thick

User benefits:

- **Controlled safety:** limits leg movement for effective containment.
- **Guaranteed comfort:** anti-friction padding to avoid skin irritation.

- **Personalized adjustment:** six eyelets per side for an optimal fit.
- **High-quality materials:** breathable cotton for prolonged use.
- **Advanced safety option:** compatible with buttons and magnetic key.

FAQ:

1. **Why is the spreader strap used?**
It is used to contain the legs in a spread position, limiting muscle movement.
2. **Who can use this fixation belt?**
Only specialized healthcare personnel can apply and remove the device.
3. **Are safety buttons included?**
No, the safety buttons and the magnetic key are sold separately.
4. **Is the material safe for the skin?**
Yes, it is made of 100% cotton with anti-friction padding for maximum comfort.
5. **Is the length adjustable?**
Yes, it has six adjustment eyelets per side and a total length of 260 cm.

Mini purchase guide:

If you need a **safe and adjustable spreader strap for healthcare environments**, the **H983_03** model is the ideal choice. Designed to **contain the legs in a spread position**, it offers **strength, comfort, and a personalized adjustment system**. Available with advanced safety closure (buttons and key sold separately).

Spreader strap for bed, leg restraint system, soft restraint devices, restraint belts for hospitals, movement limiter for patients, ankle safety strap, restraint systems with magnetic closure.

**Image purely indicative.*

**Safety buttons and safety key are available for purchase separately as options; the safety buttons can only be opened using the safety key.*

INFORMATION

