



Restraint belt for wheelchair h983_13

Belt with safety closure. Color: cream. Total length: 130 cm. Adjustable flap circumference: from 70 to 120 cm.

The **H983_13 wheelchair restraint belt** is a medical aid designed to **ensure safety and stability for patients sitting in a wheelchair**. Ideal for non-autonomous patients or those with reduced mobility, it helps **prevent slipping, accidental falls and uncontrolled movements** during their stay in a wheelchair or assistance chairs.

The device is made of **soft, resistant and breathable materials**, designed to guarantee **comfort even during prolonged use**. The **internal padding** is in polyurethane covered in MTP velvet, while the exterior and the straps (webbing) are made of **100% cotton**, gentle on the skin. The **eyelets** for the closure are in **Nickel/Brass alloy**, compatible with **magnetic safety buttons** (not included, see options).

For **advanced protection**, it is possible to purchase **safety buttons and the magnetic key separately**, which prevent unauthorized opening. The system can **only be unlocked using the dedicated key**.

To be used exclusively by **specialized healthcare personnel**.

The belt has a **total length of 130 cm**, a **height of approximately 15 cm** and an **adjustable circumference from 70 to 120 cm**, adapting to most body types. It must be worn **over clothing**, never in direct contact with the skin.

The use involves fixing it to the frame of the wheelchair backrest via the rear straps and the front closure around the patient's waist, locked with the magnetic safety system. In case of emergency, it is always necessary to have **scissors or another sharp object** available to quickly remove the device.

Protect your patients from slipping and falling with the **H983_13 restraint belt**. Order now or request a quote for multiple supplies!

Technical features:

- **Use:** patient restraint on wheelchair or seat
- **Total length:** 130 cm
- **Front height:** ~ 15 cm
- **Adjustable circumference:** from 70 to 120 cm
- **Padding material:** 100% polyamide velvet + ether-based polyurethane
- **Exterior and straps:** 100% cotton
- **Eyelets:** Nickel/Brass alloy
- **Safety closure:** with magnetic buttons and key (not included, see options)
- **Color:** cream

Optional:

- **2 safety buttons:** automatic closing buttons to be inserted into the overlapping eyelets of the safety belts
- **1 safety key:** polyamide key with 8 mm thick reinforced neodymium magnet

User benefits:

- **Prevents:** accidental slipping in the wheelchair.
- **Easy:** to apply and adjust according to the body size.
- **Soft:** breathable, ideal for prolonged use.
- **Compatible:** with magnetic safety closures.
- **Made of:** resistant and hypoallergenic materials.

FAQ:

1. What is the H983_13 belt for?

It is used to safely restrain patients sitting in wheelchairs, avoiding slipping and falls.

2. What is the size or measurement of the belt?

It is 130 cm long and adjusts in circumference from 70 to 120 cm, height approximately 15 cm.

3. Is the magnetic closure included in the kit?

No, the magnetic buttons and the key must be ordered separately (see options).

4. Can the belt be used in contact with the skin?

No, it must be worn over clothes or underwear.

5. Is it safe in case of emergency?

Yes, but it is always recommended to have scissors or a cutter available for immediate release.

Mini purchase guide:

To ensure **safety and stability** for patients in wheelchairs, it is important to choose an adjustable restraint belt, made of **breathable, comfortable and robust materials**.

The H983_13 belt offers optimal protection in hospital, nursing home (RSA) or home settings, thanks to **compatibility with magnetic safety closure systems** and **adjustable circumference**.

It is essential to respect the instructions for use, always apply the belt over clothing and have an **emergency system for quick release**.

Wheelchair restraint belt, patient restraint belt, restraint strap for disabled people, wheelchair safety belt, anti-slip aid for patients, adjustable medical belt for wheelchair, padded restraint belt for the elderly, hospital restraint belt.

**Image purely indicative.*

**Safety buttons and safety key can be purchased separately as options; the safety buttons can be opened exclusively using the safety key.*

INFORMATION

