



Electric rehabilitation pedal exerciser h1017_05

Double adjustable electric pedal exerciser for arms and legs. LCD display. Assembled product dimensions: approx. W 72 x D 63.5 x H from 95 to 104.8 cm.

The **electric rehabilitation pedal exerciser H1017_05** is a professional tool designed to support **motor recovery of arms and legs** in a safe, progressive, and customizable way. Thanks to the **dual pedaling system**, it allows you to perform **passive or active** exercises, improving **joint flexibility, blood circulation, and muscle mobility**.

It is ideal for patients in the **physical rehabilitation** phase, elderly people, or those wishing to perform light cardiovascular activities while **seated**. The reinforced structure in **cast steel, PVC, and plastic** ensures stability and durability, while the **two 70W motors** guarantee smooth and adjustable movement.

Thanks to **integrated technology**, the pedal exerciser is equipped with a **training computer with LCD display**, which shows **speed, time, distance, calories burned, and heart rate** in real time, helping the user monitor their progress. It features **5 preset programs** and a **reverse function** to vary the training type.

It is fully **adjustable in five positions**, with a **pneumatic tilt system** that allows for personalization based on the user's body size. The **pedals with adjustable straps** and the **touch-sensitive grips** (removable) promote comfortable and stable exercise. The presence of **fixing plates with stabilizers** ensures maximum safety on any type of surface.

Usable both **at home** and **in the medical field**, the h1017_05 pedal exerciser is a **versatile and compact solution**, designed to be easily moved, used anywhere, and stored without taking up space.

Buy the **h1017_05 electric pedal exerciser** today and improve motor well-being safely, even from home! Contact us for supplies to medical or rehabilitation centers.

Technical characteristics:

- **Type:** double electric pedal exerciser (arms and legs)
- **Materials:** reinforced steel, PVC, and plastic
- **Color:** gray and green
- **Motors:** 2 motors – 220V, 50/60 Hz, 70W
- **RPM:** 28/76 rpm
- **Assembled product dimensions:** approx. W 72 x D 63.5 x H from 95 to 104.8 cm
- **Net weight:** 16 kg
- **Maximum capacity:** 120 kg (seated exercise, never standing)
- **Display:** LCD with LED indicators (time, speed, distance, calories, heart rate)
- **Programs:** 5 pre-installed, with reverse function
- **Tilt adjustment:** 5 positions with pneumatic hardware
- **Grips:** removable, touch-sensitive
- **Pedals:** with adjustable straps
- **Fixing plates with stabilizers:** included
- **Timer and speed:** adjustable
- **Power cord length:** 1.5 m
- **Batteries:** not included

User benefits:

- **Allows:** training of arms and legs simultaneously or separately.
- **Suitable:** for passive or active rehabilitation.
- **Promotes:** the recovery of joint flexibility.
- **Usable:** comfortably while seated.
- **Ideal:** for the elderly, people in physiotherapy, or with reduced mobility.

FAQ:

1. **Can I use it to train only the legs or only the arms?**
Yes, the pedal exerciser is designed to work independently for arms and legs.
2. **Is it suitable for elderly people or those in a post-operative phase?**
Yes, it is ideal for light physical activity and rehabilitation, even in a medical setting.
3. **What color is the pedal exerciser?**
The pedal exerciser is gray and green.
4. **Can it be used outdoors?**
Yes, it is suitable for both indoor and outdoor use, as long as it is protected from rain and sources of humidity.
5. **Can I stand up during use?**
No, seated use is mandatory. It is not suitable for standing use.

Mini buying guide:

If you are looking for equipment for **motor rehabilitation** or to keep circulation active while **seated**, the H1017_05 pedal exerciser is an excellent choice.

Always check the available **dimensions** and ensure the patient has a safe space around them.

This model is suitable for users up to **120 kg** and is recommended for use on stable surfaces, thanks to the **non-slip stoppers**.

The **integrated programs**, the **reverse function**, and the **possibility to adjust speed and tilt** make this device perfect for both maintenance exercises and **controlled rehabilitation** paths, from home or in medical centers.

Electric rehabilitation pedal exerciser, motorized arm and leg pedal exerciser, physiotherapy pedal exerciser, electric pedal exerciser with LCD display, pedal exerciser for the elderly, motor rehabilitation tool, passive pedal exerciser for home.

**Image purely indicative.*

INFORMATION

