



Professional trolley for sterile containers h410_37

Trolley for sterile baskets with 2 vertical sections, 9 baskets each.

Transport trolley for sterile containers made entirely of anodized **AlMg3 light alloy sheet** (10 microns). **2-leaf door with 270° opening**, stop magnets in open position and vertical bolt closure with top and bottom latching. Both the central door joint and the **4 opening sides** ensure **optimal insulation of the content** during transport in non-sterile environments. **260° door locking system** for washing in sterilization tunnels. **2 external push handles**, one on each short side. **Sturdy anti-torsion base frame** in **stainless steel**. Coupling / decoupling system to the platform trolley made of **stainless steel**. Perimeter base bumper in non-marking grey rubber. **2 fixed wheels and 2 swivel wheels** diameter 160 mm in Superelastic rubber with watertight bearings and stainless steel brackets. Central brake. These cabinets are suitable for truck transport and **washing / disinfection in tunnels at 80°C**. The container is equipped with two internal sections with a structure made of anodized light alloy tubes (10 microns) and **AISI 304 stainless steel** profiles. **Shelf/basket guides** made of stainless steel. Light alloy base with **4 wheels dia. 50 mm** in fiberglass nylon for sliding in the stainless steel guides. Maximum capacity of each vertical section **9 ISO baskets 400x600x100 mm h**.

Technical features:

- Trolley for transporting shuttles for sterile containers
- Dimensions: 1230x710xh.1555 mm
- 2 shuttles for transporting ISO NORM containers 451x600x1180 mm 9 ISO NORM containers 600x400xh.100 mm

Optional:


- Transport trolley 742x510xh950 mm for 1 ISO NORM 1-section shuttle
- Unloading platform for 2 ISO NORM 1-section shuttles 1230x710xh260 mm

*image for illustrative purposes only


INFORMATION

- Length in millimetres 1230.000000
- Depth in millimetres 710.000000
- Height in millimetres 1555.000000

Professional trolley for sterile containers h410_37



HOLITY.COM



HOLITY.COM



HOLITY.COM



HOLITY.COM